

**Safe Touch® Advanced Guard Nitrile Gloves
Item #SFTGN1138X**

Item Code & Size	<table border="1"> <thead> <tr> <th>Item Code</th> <th>Size</th> </tr> </thead> <tbody> <tr> <td>SFTGN1138A</td> <td>Extra Small</td> </tr> <tr> <td>SFTGN1138B</td> <td>Small</td> </tr> <tr> <td>SFTGN1138C</td> <td>Medium</td> </tr> <tr> <td>SFTGN1138D</td> <td>Large</td> </tr> <tr> <td>SFTGN1138E</td> <td>Extra Large</td> </tr> </tbody> </table>		Item Code	Size	SFTGN1138A	Extra Small	SFTGN1138B	Small	SFTGN1138C	Medium	SFTGN1138D	Large	SFTGN1138E	Extra Large
	Item Code	Size												
	SFTGN1138A	Extra Small												
	SFTGN1138B	Small												
	SFTGN1138C	Medium												
	SFTGN1138D	Large												
SFTGN1138E	Extra Large													
Country of Origin	Malaysia													
Length	Min 240mm													
Width	Extra Small (XS)	Min 73 mm												
	Small (S)	Min 81 mm												
	Medium (M)	Min 91 mm												
	Large (L)	Min 102 mm												
	Extra Large (XL)	Min 110 mm												
Thickness	Min 0.06 mm (Palm) Min 0.10 mm (Finger)													
Tensile Strength	Before aging ≥ 14 MPa After aging ≥ 14 MPa													
Elongation %	Before aging ≥ 500% After aging ≥ 400%													
Colour	Black													
Texture	Finger Textured													
Sterility and usage	Non-sterile, Disposable, Single-Use only													
Powder Residue	≤ 2mg per glove													
Packaging	100 gloves per box / 10 boxes per case / 1000 gloves per case													
Storage	Store under cool, clean, and dry conditions. Avoid excessive heat (over 40°C or 104°F). Shield open box from direct sun and fluorescent lighting													
Shelf-life	5 years													

This glove meets ASTM Standard D 6319
--

This technical data sheet issued by A.R. Medicom Inc (Asia) Ltd